

DDM

DESSILLONS & DUTRILLAUX - Since 1949

TECHNICAL CHARACTERISTICS

CT 059

Pressure infusion cuff reusable Grey
Fuse

Index 5

MADE IN FRANCE **Regulatory information**

CE Marking	Index of classification	Notified body	PHTALATES	LATEX
CE 0459	Ila	G-MED		

1 – Device identification

The Grey Fuse pressure infusion cuffs are made from one set :

- robust design
- High quality technical fabric
- cleaning and decontamination facilities
- strap for hanging the pressure infusion cuff (on M25000 only)
- hook to keep the infusion bag on M25000 due of weight of fluid bags
- The transparent monofilament view window tear-proof to view the liquid

**2 - References**

Reference	Designation	Size	Dimensions in mm	
			Lenght	Width
M20500	Monobloc pressure infusion cuff	500 cc	250	162
M21000	Monobloc pressure infusion cuff	1 000 cc	365	170
M25000	Monobloc pressure infusion cuff	3000 to 5 000 cc	561	350

Dessillons & Dutrillaux

Zone Industrielle La Tuque
Tél : +33 5 53 48 30 66
technique@ddmedical.fr

47240 CASTELCULIER
Fax : +33 5 53 47 24 64
www.ddmedical.fr

Date of first EC marking certificate :
2004


Page 1 / 4

3 – Use of the device - Operating instructions

Visual inspection : before use the medical device.


- On a horizontal plane, positioned under the transparent fabric, the pocket to press, output tubes oriented to pressure gauge.
- Insert the fastening system of pressure cuff in the opening of the pocket to be irrigated, before slipping into the opening made on the pressure cuff.
- Using the bulb, breathe the pressure required to maintain the cuff around the bag to irrigate, suspend all on a drip stand. Inflate the pressure cuff until the desired pressure (max : 300 mmHg). In case of exceeding pressure greater than 300 mmHg, the red zone will appear and the pressure will be regulated automatically.
- As to the flow of product irrigated, the pressure displayed by the pressure gauge falls gradually : it is necessary to restore order to ensure complete emptying of the liquid.
- The decompression of the cuff is facilitated by the clip connectors.

4 - Identification and characteristics of different subsets



Identifier	Reference	Denomination	Characteristics
I	M20511	Pressure cuff only 500 cc	PU coated fabric+ poyamide monofilament
I	M20011	Pressure cuff only 1000 cc	PU coated fabric+ poyamide monofilament
I	M25011	Pressure cuff only 3000 - 5000 cc	PU coated fabric+ poyamide monofilament
O	A10808	Pressure gauge	ABS PC +glass PC + TPU
C	A10259	Valve bulb	ABS
A-B-C	A10226	Complete bulb	PVC phtalates free
F	A11012	Female coupler online fluted 4.8 mm	ACETAL
G	A11008	Male coupler online fluted 4.8 mm	ACETAL EPDM
E	A10130	Connectors fastener (A11012 + A11008)	TPU

Identification of the pressure gauge components :



DDM

DESSILLONS & DUTRILLAUX - Since 1949

TECHNICAL CHARACTERISTICS

CT 059

Pressure infusion cuff reusable Grey
Fuse

Index 5

5 – Cleaning and disinfection :

Surface cleaning	Disinfection with decontamination solution
Complete device	Complete device

6 - Maintenance :

Visual inspection / watertight : before each use.
Calibration test : every 12 months.

7 - Packaging

Reference	Packaging units
M20500	1
M21000	1
M25000	1

8 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10°C à 40 ° C	30 à 40 %	500 à 1060 hPa

9 - Warranty :

This warranty provides assurance for customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 2 (two) years from the date of purchase.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Desillons & Dutrillaux

Zone Industrielle La Tuque
Tél : +33 5 53 48 30 66
technique@ddmedical.fr

47240 CASTELCULIER
Fax : +33 5 53 47 24 64
www.ddmedical.fr

Date of first EC marking certificate :
2004

Page 4 / 4