


MADE IN FRANCE 

1 – Device identification

CE Marking	Index of classification	PHTALATES LATEX
CE	I	



Cuff are used for blood pressure with sphygmomanometer. Arm cuff are made of polyurethane coated fabric and they are immersible.

Flexible and robust design, it's possible to clean in the washing machine and decontaminate them by immersion with the obligation to close the connection tubes.

- Reliability : an Index and a range "Min-Max" can quickly validate if the product is adapted to the size of the patient.

- Security : to fight against nosocomial infections, this fabric underwent an antimicrobial treatment.

2 - References

Reference	Designation	Tube	Circumference in cm		Dimensions in mm					Tube color
			Min	Max	A	B	C	D	E	
BM101	New born	1	6	12	170	40	95	47,5	43	WHITE
BM102	Child	1	11	22	338	80	180	90	85	PINK
BM103	Small adult	1	16	28	423	110	230	115	120	PARMA
BM104	Adult	1	26	33	516	129	256	129	70	BLUE SKY
BM105	Adult large	1	33	41	616,8	150	330	165	80	TURQUOISE
CM104	Lower limb cuff adult	1	39	55	816,5	190	440	220	160	VIOLET
BM201	New born	2	6	12	170	40	95	47,5	43	WHITE
BM202	Child	2	11	22	338	80	180	90	85	PINK
BM203	Small adult	2	16	28	423	110	230	115	120	PARMA
BM204	Adult	2	26	33	516	129	256	129	70	BLUE SKY
BM205	Adult large	2	33	41	616,8	150	330	165	80	TURQUOISE
CM204	Lower limb cuff adult	2	39	55	816,5	190	440	220	160	PURPLE

Dessillons & Dutrillaux

Zone Industrielle La Tuque

Tél : +33 5 53 48 30 66

technique@ddmedical.fr


47240 CASTELCUJER

Fax : +33 5 53 47 24 64

www.ddmedical.fr


Date du premier certificat de
marquage CE : 2004












Page 1 / 5



3 – Technicals characteristics

- Latex free PU coated fabric - high-frequency welding
- Loop Astrakan high-frequency welding



A		Indication of the position of the index
B		Product Brand
C		Name and product number
D		Original product marking
E		Range of product use
F		Date of manufacture
G		See operating instructions
H		Latex free
I		Manufacturing Lot Number
J		Validity range of product use
K		Indication for placement of the artery

4 – Operating instructions

- Select a cuff adapted to the morphology of the patient. The size and range of use are indicated on the cuff.
- Position the center of the bladder pocket represented by a heart on the axis of the brachial artery.
- Roll the cuff around the limb.
- Once the cuff positioned and ready to use, the INDEX must be between the Min and Max (axis of the heart). If INDEX exceeds the Max index located at the axis of the heart, replace the cuff by the model of the larger size. If INDEX exceeds Min index located at the tip of the range, replace the cuff by the model of smaller size.

- When the cuff and INDEX are positioned correctly the procedure for measuring blood pressure can begin.
- Correctly roll the cuff around the limb and ensure its continuation with the fastening system provided.
- Place the tips of the lyre of the stethoscope in your ears, the microphone in the listening position, just below the cuff, on the passage of the artery.
- Ensure that the vent screw of the pear is in the closed position 30 and to inflate or 40 mmHg above the estimated pressure. Blood no longer circulates at the microphone, no sound is heard with a stethoscope.
- Loosen the vent screw the bulb to slowly reduce the inflation pressure and allow blood to flow again, the first beats are perceptible stethoscope.
- Read this pressure gauge reading now, this value corresponds to the maximum pressure or systolic pressure.
- Continue decompression beats levied stethoscope amplify then fade to become inaudible.
- Read this moment the pressure indicated, this reading corresponds to the minimum pressure or diastolic pressure.
- Unscrew the knob further relief to completely purge the remaining air in the bladder.
- Remove the cuff.

5 - Maintenance :

N A

6 - Cleaning :

1) Dry cleaning :

- A) Spray a detergent, disinfectant, bactericidal foam on the arm cuff, taking care to distribute the product well on the surface.
- B) Leave for 15 minutes.
- C) Do not rinse.
- D) Product recommended: Surface'Safe / Anios.

	INSTRUCTIONS FOR USE		INS 104
	Immersible arm-cuff "Easy Cuff"		Index 4

2) Cleaning by immersion :

- A) Seal the tube or tubes by means of the closure system provided :
- automatic valve : ref A11300 ; cuff's side + A11308
 - Connector with cap: A11360 ; cuff's side + A11302
 - Arm cuff with dual output, connect the tubes together.
- B) Clean the arm cuff with water or with washing machine. Washing machine at 40 ° (60 ° being the maximum), Without a spin cycle. The number of washing machine and the washing temperature decreases the life of the cuff.
- C) Immerse the tube and arm cuff assembly in a decontaminating solution.
- D) Immersion time : 15 minutes.
- E) Rinse thoroughly with water, renewing the dipping bath after each use. Product recommended: Aniosyme DD1 manufactured by Anios.

7 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

8 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 1 (one) year from the date of purchase.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Dessillons & Dutrillaux Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr		47240 CASTELCUJER Fax : +33 5 53 47 24 64 www.ddmedical.fr	Date du premier certificat de marquage CE : 2004	Page 5 / 5
-----------------------------------------------------------------------------------------------------------------------	--	------------------------------------------------------------------	-----------------------------------------------------	------------