

# DO MORE FOR LESS.

## SenSmart® 8006 Hybrid Sensor



### True Innovation Improves Performance and Reduces Cost

Nonin's SenSmart® 8006 Hybrid Sensor does both. Clinically proven dual emitters and absolute  $rSO_2$  measurements improve accuracy. A reusable sensor cuts annual spend by up to 50% while reducing waste.\* And a disposable patient interface helps prevent cross-contamination. Introducing SenSmart 8006, the next generation  $rSO_2$  sensor.

**Innovation without compromise.**

\*Based on Nonin economic value calculations



# Introducing the Next Generation rSO<sub>2</sub> Sensor

## Advance Your Clinical and Cost Effectiveness



### Innovative

Dual emitter, 16 light path rSO<sub>2</sub> sensor reduces signal contamination for absolute accuracy



### Disposable

Single-use patient interface helps prevent cross-contamination; Nonin's reliable adhesive ensures sensor stays in place



### Affordable

Up to 50%\* average reduction in annual sensor spend



### Durable

Designed and tested to last for a minimum of 300 uses\*\*; yields significant cost savings and reduces waste

\*Based on Nonin economic value calculations  
\*\*Data on file at Nonin, sensor has 6 month warranty

## Reliable Measurements Clinicians Trust

The SenSmart 8006 Hybrid Sensor works with the SenSmart Universal Oximetry System.

Proven SenSmart rSO<sub>2</sub> technology provides reliable and consistent measurements

Small footprint provides room for other monitoring sensors

Versatile design allows use on any patient > 40kg

rSO<sub>2</sub> and SpO<sub>2</sub> values appear on the same screen for real-time decision making

SenSmart system allows monitoring of up to six sites on a patient

Indicator light on sensor confirms proper connection



### Nonin Medical, Inc.

13700 1st Avenue North  
Plymouth, MN • 55441-5443 • U.S.A.

Tel: +1.763.553.9968 1.800.356.8874

Fax: +1.763.577.5521

Email: info@nonin.com

### Nonin Medical B.V.

Doctor Paul Janssenweg 150  
5026 RH Tilburg • Netherlands

Tel: +31 (0)13-79 99 040

Fax: +1 (0)13-79 99 042

Email: infointl@nonin.com



nonin.com



©2020 Nonin Medical, Inc. All trademarks are the property of Nonin Medical, Inc. unless otherwise noted.  
M-18-010-05

CE 0123



UL 00001-1  
2020M

Rx Only