


**HFOT2000**

**HFOT2000**  
operation position sketch



**Main Technical Data**

HFOT2000 Manual Hydraulic Operating table	
Tabletop Length	2085mm
Tabletop Width	550mm
Elevation (down/up)	520/820mm
Trendelenburg	≥25°
Reverse Trendelenburg	≥25°
Lateral Tilt(left/right)	15°/15°
Head Rest(up/down)	40°/90°
Footplate(up, down, outward)	15°/90°/90°
Back plate(up/down)	80°/20°
Kidney bridge elevation	110mm
Power	220V 50Hz 1.0kW

