




HF EOT99C

HF EOT99C
operation position sketch



Main Technical Data

HF EOT99C	
Tabletop Length	1970mm
Tabletop Width	500mm
Elevation (down/up)	690/1100mm
Trendelenburg/Reverse Trendelenburg	30°/30°
Lateral Tilt(left/right)	20°/20°
Head Rest(up/down)	40°/90°
Foot Plate(up/down& Horizontal Positioning)	15°/90°/90°
Back Plate(up/down)	90°/40°
Longitudinal Sliding	350mm
Power	220V 50Hz 1.0kW
Kidney bridge elevation	110mm

