

HFMED[®]

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

HFEOT99S Electro-Hydraulic Operating Table

Developments made in the conjunction with doctors; this is one of the basic principles which has made the world market leader for operating tables. Close cooperation and years of experience have resulted in tried and tested products which facilitate the work in the operating room and help to increase the safety and comfort of the patient.

Best example: this table can maintain stability in any position even with extreme loads. The modular structure enables it to be adapted to match the size of the patient, and in no time at all it can be equipped for any surgical discipline. A longitudinal shift of up to 350mm makes the operating table even more versatile.

Your Ideal Operating
Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery



HFMED[®]

SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD

Address: NO.885, Qiuxing Road, Nicheng Town, Pudong New Area, Shanghai, China.
Tel: 13764286262



HFMED

YOUR IDEAL OPERATING TABLE SUPPLIER
MEET THE VARIOUS REQUIREMENTS OF GENERAL SURGERY

HFEOT99S

Electro-Hydraulic Operating Table



SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

HFEOT99S Electro-Hydraulic Operating Table



Product Description

Engineered for increasingly complex surgical scenes, HFEOT99S comes with intelligent automatic control system satisfying various demands of general surgeries. Extra wide tabletop, extreme loading system, long horizontal sliding that can be suitable for both X-ray and C-arm. HFEOT99S is an ideal electro-hydraulic operating table.

- Longitudinal slide up to 350mm for free access to C-arm
- Extreme weight load capacity of 300kg with highest safety and stability
- Memory pad with artistic, waterproof design.
- Standard override panel control makes the control safer
- "EASY CLICK" leg plate control by electro-hydraulic and manual
- Battery power for 50-80 operations



PARAMETER HFEOT99S Electro-Hydraulic Operating Table

Movements		Model: HFEOT99S
Tabletop: Trendelenburg	Reverse Trendelenburg	≥30°
	Trendelenburg	≥30°
Lateral Tilt (Left And Right)		≥20°
		≥60°
Head plate	Up	≥90°
	Down	≥90°
Back plate	Up	≥40°
	Down	≥15°
Leg plate	Up	≥90°
	Down	≥90°
Kidney bridge elevation		110mm
	Horizontal sliding	≥350mm

Details



1. Maximum 350mm motorized longitudinal shift for flexible use of the C-arm.
2. "EASY CLICK" for simple and quick module changes.
3. Leg plate imported gas spring for easy controls.
4. One button reposition function, designed and made specially for X-ray and C-arm use.
5. Double control system.
6. Built-in kidney bridge.
7. Two-joint head plate.



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

Model	Site	Length	Width	Tabletop lowest height	Elevation range
HFEOT99S		≥2170mm	550mm	≤695mm	≥410mm

STANDARD ACCESSORIES

Anaesthesia screen	1 piece	Body support	1 pair	Leg support	1 pair	Hand remote	1 piece
Kidney bridge handle	1 piece	Arm support	1 pair	Foot support	1 pair	Matress	1 set
Shoulder support	1 pair	Long fixing clamp	2 pieces	Fixing clamp	7 pieces	Power Line	1 piece

