


HFMH3008AB Manual Hydraulic Operating table

HFMH3008AB operation position sketch



Main Technical Data

HFMH3008AB Manual Hydraulic Operating table	
Length and width for delivery	2100/900mm
Elevation(down/up)	780/1020mm
Trendelenburg/Reverse Trendelenburg	22/22°
Lateral (left/right)	18/18°
Head rest(up/down)	45/70°
Back plate(down)	75/115°
Leg plate(down/outward)	15/90°/0°
Maximum load capacity	180kg

