

PRODUCT MODEL



HF-LED SERIES TECHNICAL DATA

Type	Illuminance	Colour temperature	Colour reduction index(Ra)	Total irradiance	Size of light field	Service life of illuminant	Power supply voltage	Best height of installation	Total power consumption	Total LED bulb quantity	Green ENDO	Dimming Type
HF-L3+4C LED	16000lx 14000lx	3500±300k 4000±300k 4500±300k 5000±300k	Ra:96 R9:97 R13:98	380 380	120-250mm	50000h	AC110-240V 50-60Hz	2900mm	90W	187pcs	YES	Analog dimming
HF-L4+4C LED	16000lx 16000lx	3500±300k 4000±300k 4500±300k 5000±300k	Ra:96 R9:97 R13:98	380 380	120-250mm	50000h	AC110-240V 50-60Hz	2900mm	100W	224pcs	YES	Analog dimming
HF-L3+3C LED	14000lx 14000lx	3500±300k 4000±300k 4500±300k 5000±300k	Ra:96 R9:97 R13:98	380 380	120-250mm	50000h	AC110-240V 50-60Hz	2900mm	80W	150pcs	NO	Analog dimming
HF-L3C LED	14000lx	3500±300k 4000±300k 4500±300k 5000±300k	Ra:96 R9:97 R13:98	380	120-250mm	50000h	AC110-240V 50-60Hz	2900mm	40W	75pcs	NO	Analog dimming
HF-L4C LED	16000lx	3500±300k 4000±300k 4500±300k 5000±300k	Ra:96 R9:97 R13:98	380	120-250mm	50000h	AC110-240V 50-60Hz	2900mm	50W	112pcs	YES	Analog dimming